




Ethicon Products


W O R L D W I D E


MERSILENE

**Braided and Monofilament
Polyester Suture**


Braided Polyester Suture



Needle length and description	Length	Colour	Gauge						
▼ 3/8 circle Reverse Cutting			3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	5 7
 FS-1	24mm	45cm	white	R647H					


Needle length and description	Length	Colour	Gauge						
▼ 1/4 circle SPATULATED			5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
 S-14	8mm	45cm	white	W843	W832				


▼ 1/2 circle SPATULATED			5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
 S-22	8mm	45cm	white	W833					

MERSILENE Monofilament Polyester Suture

Needle length and description	Length	Colour	Gauge						
3/8 circle CS ULTIMA spatula			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5
 CS140-6	6.5mm	30cm	green	U7075					

▼ 3/8 circle Advanced MICRO-POINT spatula			11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 TG140-6	6.5mm	30cm	green		W1707				
 TG140-6	6.5mm	30cm	green	W1785	W1775				

▼ 1/2 circle Advanced MICRO-POINT spatula			11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 TG160-6	5.5mm	30cm	green	W1758	W1755				

▼ Compound Curve Advanced MICRO-POINT spatula			11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 TG6-C	5mm	30cm	green	W7795					

Monofilament Polyester Suture



Needle length and description		Length	Colour	Gauge	
CATARACT PROCEDURE PACK				10-0 0.2	
ADVANCED MICRO-POINT spatula 3/8 circle		TG140-6  6.5mm	25cm	green	W7720
Each pack also contains one 4-0 black MERSILK braided silk suture, 30cm, attached to a 16mm 3/8 circle cutting needle.					

Ligatures and Non-needed Sutures

TROCHANTER SUTURE

Gauge Metric	white 1.5m length 12 sachets per box
9	W394

MERSILENE Tape

Needle length and description		Length & Colour	Width of tape	
48mm				
1/2 circle round bodied heavy double	● 	40cm white	5mm	RS22
65mm				
BLUNT POINT 1/2 circle round bodied heavy double	● 	30cm white	5mm	RS21



Capitalised product names are trademarks of ETHICON.