












Ethicon Products









W O R L D W I D E

**SURGICAL
STAINLESS STEEL
WIRE**

Monofilament and Multifilament

Needle length and description		Length	Colour	Gauge						
☉ 1/2 circle Taper Point				4-0 2	3-0 3	2-0 3	0 3.5	1 4	2 5	4 6
 MH-1	31mm	75cm	multi	W310	W311					
▼ 1/2 circle TAPERCUT				2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9
 V-37 heavy	40mm	4x45cm	mono			M660G				
STRAIGHT PACK										
 V-40 heavy	48mm	4x45cm	mono					M651G	M650G	
STRAIGHT PACK										
▼ 1/2 circle Reverse Cutting				2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9
 CPX	48mm	4x45cm	mono							M624G
STRAIGHT PACK										
▼ 1/2 circle Reverse Cutting Rotating Needle				2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9
 CPXX	55mm	4x45cm	mono							M635G
STRAIGHT PACK										
 CPXX	55mm	75cm	mono							W905
STRAIGHT PACK										
 CPXX	55mm	75cm	mono							W995

Needle length and description	Length	Colour	Gauge						
▼ 1/2 circle Reverse Cutting Rotating Needle			2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9
CPXX  STRAIGHT PACK	55mm	2x75cm	mono						W902
▲ 1/2 circle TROCAR POINT			2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9
TR-55 	55mm	75cm	mono						W945

Needle length and description	Length					
		5-0 1	4-0 1.5	3-0 2	2-0 3	1 4
PULL-OUT TENDON SUTURES						
<p>A complete suture set for tendon repair of the hand. Each set contains a double-armed suture, a wire loop swaged to a needle and a button. This product has been designed for the Bunnell technique.</p> <p>26mm 3/8 circle TAPERCUT </p> <p>51mm straight cutting double </p>	<p>20cm loop Mono</p> <p>35cm Mono</p>					D9805
PULL-OUT TENDON SUTURES						
<p>50mm straight Taper Point </p> <p>22mm 1/2 circle reverse cutting </p>	<p>30cm Multi with barb</p>					FSN13
TEMPORARY CARDIAC PACING WIRE <i>insulated with blue polyethylene.</i>						
<p>Incorporating a 60mm straight needle which has a 'breakaway' feature enabling the remaining butt to be used as a connecting plug to the Pacing Unit.</p>						
						Diameter
						0.3mm
<p>17mm TAPERCUT 1/2 circle </p> <p>60mm straight cutting </p>	<p>60cm Multi</p>					FEP13E
<p>26mm TAPERCUT 1/2 circle </p> <p>60mm straight cutting </p>	<p>60cm Multi</p>					FEP15E