




Ethicon Products

W O R L D W I D E

BLAKE
Silicon Drains



Virtually pain-free on removal

Uniquely designed to be patient friendly, **BLAKE* Silicone Drains** are an innovative surgical drain with a unique 4-channel design to provide greater tissue-to-drain contact and multiple drainage routes. Constructed of silicone with a solid centre core, **BLAKE Drains** are smooth, small and flexible, allowing easy insertion and positioning. With no holes for the tissue to invaginate into, **BLAKE Drains** will minimise pain on removal, sliding easily and smoothly past any tissue ingrowth.


For maximum patient comfort and ease of use, choose **BLAKE Drains**.

ETHICON BLAKE Silicone drains

BLAKE Silicone drains


Round BLAKE Hubless Silicone Drains

Code	Description	Packaging
2226	10Fr. Round Hubless	10 per box
2227	10Fr. Round Hubless with 1/8 inch Trocar	10 per box
2228	15Fr. Round Hubless	10 per box
2229	15Fr. Round Hubless with 3/16 inch Trocar	10 per box
2233	15Fr. Round Hubless with 3/16 inch Bendable Trocar	10 per box
2231	19Fr. Round Hubless with 1/4 inch Trocar	10 per box
2232	19Fr. Round Hubless with 1/4 inch Bendable Trocar	10 per box
2230	19Fr. Round Hubless	10 per box
2234	24Fr. Round Hubless	10 per box



Flat BLAKE Silicone Drains-Full Fluted


2214	10mm Flat	10 per box
2215	10mm Flat with 3/16 inch Trocar	10 per box



J-VAC Reservoirs


J-VAC Bulb Suction Reservoir

2160	100cc Reservoir	10 per box
------	-----------------	------------



J-VAC Reservoir

2161	150mL Reservoir	10 per box
2162	300mL Reservoir	10 per box
2163	450mL Reservoir	10 per box



BLAKE Cardio Connectors

BCC1 BLAKE Cardio Connector 1:1


3/8" - 1/2" x 3/16"	(10 - 13mm x 5mm)	20 per box
---------------------	-------------------	------------

BCC2 BLAKE Cardio Connector 2:1

3/8" - 1/2" x 3/16" x 3/16"	(10 - 13mm x 5mm x 5mm)	20 per box
-----------------------------	-------------------------	------------

BCC3 BLAKE Cardio Connector 3:1

3/8" - 1/2" x 3/16" x 3/16" x 3/16"	(10 - 13mm x 5mm x 5mm x 5mm)	20 per box
-------------------------------------	-------------------------------	------------





Capitalised product names are trademarks of ETHICON.